

# CHEMTRIX

## Opportunities for the Development of Sustainable Production Processes

Dr Charlotte Wiles

RSC Waste Not Want Not Symposium, Budapest, Hungary


18-19<sup>th</sup> June 2014

Strategic Partner of



# Chemtrix Company History

**CHEMTRIX**



More than 13 years experience in Flow Chemistry

Cooperation with DSM since 2012 together offering from concept to delivery solutions

# Innovative Technology: Scalability & System Flexibility

**CHEMTRIX**



**Labtrix<sup>®</sup>** ( $\mu\text{g}$  to  $\text{mg}$ 's)  
-20 to 195 °C

**DISCOVERY**

- Rapid reactions
- Efficient evaluation
- mg consumption
- Parameter accuracy
- New chemical entities



**KiloFlow<sup>®</sup>** ( $\text{g}$  to  $\text{kg}$ 's)  
-15 to 150 °C

**DEVELOPMENT**

- Rapid up-scaling
- Process validation
- kg Production in a lab
- New process windows
- Flexible production



**Plantrix<sup>®</sup>** ( $\text{kg}$  to  $\text{ton}$ 's)  
-30 to 200 °C

**PRODUCTION**

- Facile up-scaling
- Forbidden chemistry
- Safe production
- High quality products
- Cost effective

Customised solutions are also delivered in partnership with our Customers



# Conventional Synthetic Methodology: Challenges and Limitations

CHEMTRIX

- If we look to how synthetic chemistry has been taught and performed, little has changed over the past century, with all chemists being familiar with **standard glassware** and equipment

## Batch Reactions:

- In batch reactions parameters such as time, temperature, stoichiometry, order of addition and solvent are investigated with the aim of increasing yield and product purity
- If more product is required then a **larger vessel** is normally employed



- Changes in **surface to volume ratio** mean that differences in thermal and mass transfer occur and reactions often need to be re-optimised

# Conventional Synthetic Methodology: Challenges and Limitations (2)

**CHEMTRIX**

'Batch' has been used for 100's years, so it must be good for something;

## Advantages:

- Flexible
- Known scale-up path
- Good all round approach – liquids, solids, gases tolerated

It is however known to have its limitations;

- Mixing
- Mass transfer
- Heat transfer
- Pressure limits
- Temperature limits
- Material incompatibility
- Cost – Maintenance & labour intensive




'Batch has its issues, but time has taught ways around the problems'

# Conventional Synthetic Methodology: Challenges and Limitations (3)

Whilst their flexibility means that as Chemists we persevere with batch vessels;

- Broad temperature, time & concentration profiles – leads to difficulties in reaction control

Here is illustrated a reaction with a competing by-product that consumes the target product



As the size increases, control is more problematic

- By obtaining a narrower temperature, time and/or concentration profile, it is possible in flow to prepare the target product without competing by-product formation


# Fundamentals of Flow Chemistry:

## How are Flow Reactions Performed?

CHEMTRIX

Solutions of reagents are pumped into a reactor, where they are;

- Heated or cooled ahead of mixing
- Reacted for specified period of time - based on volume of reactor & flow rate





- Quenched *in-situ* (where needed)
  - To stabilise product & prevent decomposition
- Collected for batch isolation & purification (where needed)

# Fundamentals of Flow Chemistry: How are Flow Reactions Performed? (2)

CHEMTRIX

- Reactor volumes can range from  $\mu\text{l}$ 's to  $\text{ml}$ 's, consequently;
  - ✓ Small volumes of reagents & catalysts used for large number of reactions
  - ✓ Steady-state is reached with low volume consumption
  - ✓ Reaction conditions can be changed rapidly
- Molecules experience similar conditions leading to increased process stability
- Volume production is **time-resolved**, increasing the volume of product can be achieved by increasing the run time or by employing multiple reactors or by increasing reactor volume
  - ✓ Theoretically no failure to scale and no changes in the safety profile of a process

Reactors can be micro structured or simple tube reactors depending on the application






# Fundamentals of Flow Chemistry: Efficient Heat Exchange & Mixing

CHEMTRIX

## Heat Exchange:


- Due to the reduced size of a reaction channel uniform temperatures are obtained
- Reduces the risks associated with exotherms
- Increases product quality
- Allows previously uncontrollable process to be executed



## Mixing:

- Mixing occurs by diffusion only and is linked to the channel size and the reactants employed

$$t \sim \frac{d^2}{D}$$



- In a 300 µm channel this can take 10's sec
- Incorporation of static mixers can reduce the mixing time further to the ms range

# Fundamentals of Flow Chemistry: Examples of Micro Mixers

**CHEMTRIX**


100  $\mu\text{l min}^{-1}$  (left) T-mixer & (right) SOR-mixer



# Flow Technology: Drivers for Implementation



If we look to the drivers associated with the development of a process in flow;



Source: Pharmaceutical Research and Manufacturers of America.

## Lab-scale:

- Speed
- New reaction space
- Reproducibility
- Selectivity
- Flexibility

## Process R&D:

- Speed
- Safety
- Robustness

## Production:

- Speed
- Safety
- Robustness
- Cost reduction
- Quality

**As each stage has different drivers – different equipment is used**

# Flow Technology: Sectors & Users




Applications reported in;

- Pharmaceutical, fine chemical, agrochemical, specialty & commodity industries
- Research & Development
  - New chemical processes evaluated
    - Detailed investigations performed using small quantities of reagents and catalysts
  - Material production
    - Nanoparticles, colloids, pigments , polymers prepared with high specifications
- Process Development
  - Rapid translation of processes from R&D to pilot scale
  - Can stay in the research laboratory for longer to product kilograms of material
  - Can use previously inaccessible approaches for the production of materials
    - Biocatalysis, photochemistry/electrochemistry accessible at scale
- Production
  - Previously forbidden chemistries can be performed at scale with improved safety profile
  - Synthesis of fine chemicals, pigments and ionic liquids
  - Active pharmaceutical intermediate (API's) production


# Continuous Azidation using Labtrix<sup>®</sup>: Manipulation and Formation of Hazardous Materials

**CHEMTRIX**

- Employing 0.66 M (EtOH) alkyl precursor and 0.66 M (50:50 aq. EtOH) NaN<sub>3</sub>, the effect of temperature on the reaction was investigated at a residence time of 30 s



'Novel Operating Window'  
at the mg-scale



Decreasing  
material costs

Limit of conventional  
batch glassware



- Using OMs derivative, azide obtained at a throughput of 79 mg h<sup>-1</sup> @195 °C

# Continuous Nitrations using Labtrix<sup>®</sup>: Manipulation of 70 % Nitric Acid

**CHEMTRIX**

## Reaction Conditions:

- Solvent-free & mixed acid feeds
- Equal flow rates
- Reaction times = 7.5 to 60 s
- Reactor temperature = 0 to 40 °C



HNO <sub>3</sub> :H <sub>2</sub> SO <sub>4</sub>	HNO <sub>3</sub> :hexanol (eq.)	Product	By-product
1:0	3.1	✗	✗
1:0.286	2.5	✗	✓
1:0.767	1.5	Minimal	✓
1:1.130	1.25	✓	✓
1:1.726	1.0	✓	✗

✓ No by-product formation observed under optimal conditions

## Optimal Conditions:

- 60 s reaction time @ 40 °C → System throughput = 336.3 mg hr<sup>-1</sup>

# Facile Up-scaling from Labtrix<sup>®</sup> to KiloFlow<sup>®</sup>: No Re-optimisation Required

**CHEMTRIX**

KiloFlow<sup>®</sup> is a modular, **scalable** flow reactor system that can support up to **Phase 3**

- Giving you access to a **pilot plant** within a standard laboratory fume hood

KiloFlow<sup>®</sup> Product Portfolio has a;


- Large working range
  - **-15 to 150 °C**
- Flexible production volume of
  - 0.012-6.0 l h<sup>-1</sup> (up to **140 l day<sup>-1</sup>**)
- Small footprint
- Glass meso reactors
  - Low pressure drop
  - Efficient mixing
  - Excellent heat exchange ( $U \times (S/V) = 3265 \text{ kW/m}^3 \cdot \text{K}$ )
- Allows method transfer from Labtrix<sup>®</sup>



# Production using KiloFlow<sup>®</sup> Basic: Synthesis of 1-(2-Methyl-2*H*-chromen-3-yl)ethanone

**CHEMTRIX**

Finding application in the pharmaceutical, agrochemical & flavours/fragrance sectors, coumarins are of significant interest to researchers




Using the condensation of salicylaldehyde and ethyl acetoacetate in the presence of an organic base, the synthesis of 3-acetylcoumarin was optimised in Labtrix<sup>®</sup>

## Optimal Conditions:

- 60 s reaction time
- 125 C reactor temperature
- 0.4 eq. piperidine
- MeCN reaction solvent

1  $\mu$ l Reactor volume (3221) consumed 5.4 mg h<sup>-1</sup>






# Production using KiloFlow<sup>®</sup> Basic: Synthesis of 1-(2-Methyl-2*H*-chromen-3-yl)ethanone

CHEMTRIX

Translating the optimal conditions to KiloFlow<sup>®</sup> Basic (Reactor Volume = 13 ml)



Operating for 5.2 h, **369.6 g** of 1-(2-methyl-2*H*-chromen-3-yl)ethanone was obtained

- After aq. extraction (98.2 % yield)

# Production using KiloFlow<sup>®</sup> Basic: Synthesis of 1-(2-Methyl-2*H*-chromen-3-yl)ethanone

**CHEMTRIX**



**13,000 x scale-up** without;

- Parameter re-optimisation
- Change in product quality

Turn-key Flow Platform for;


- mg-scale optimisation to Kg production



# Customer Appraisal of KiloFlow<sup>®</sup>: Janssen Pharmaceutica NV

**CHEMTRIX**

## Eschweiler Clarke Reaction:



Modelling predicted the optimal to be 56 s @ 122 °C

- Outside the conditions safely attainable in batch (time, temp & press), CO<sub>2</sub> ↑



	Batch	Predicted Optimal Conditions <sup>A</sup>	KiloFlow <sup>®</sup>
N-Methyl Derivative	87.0	95	93.5
By-product	10.0	<3.5	2.10
Others	3.0	1.2	1.0-1.2

<sup>A</sup> Predicted at 56 s @ 122 °C

- Reaction translated to a 1 litre PFA reactor for 100 kg h<sup>-1</sup> production of the API intermediate
- 3 Validation batches performed, results presented to the FDA who confirmed;

**'no additional analytical PAT tools were required for production'**

# Customer Appraisal of KiloFlow<sup>®</sup>: Iolitec GmbH



**CHEMTRIX**

## Chemical Appraisal:


- Continuous flow synthesis of Ionic Liquids (IL's)
  - Uses include solvents, process chemicals, functional fluids, electrolytes, additives
- IL is a general term for a class of materials consisting of ions and being liquid below 100 °C (RTIL's) consequently there are 10<sup>18</sup> potential combinations of ions

## Iolitec have a diverse range of products (~ 300 IL's):

- Production quantities range from g to tonnes
- Rapid and scalable production techniques are required for high purity IL synthesis

## Solution: Employ KiloFlow<sup>®</sup> with integrated heat exchangers

- Higher HE capacity than competitors; **3265 kW/m<sup>3</sup>.K**



# Plantrix<sup>®</sup> made of EKasic<sup>®</sup> Silicon Carbide


## Efficient Industrial Production & Superior Chemical Flexibility

**CHEMTRIX**



- High productivity
- High chemical flexibility
- Less scale-up risk
- Handling of solids
- Increased safety
- Environmental friendly production
- Small footprint

# EKasic<sup>®</sup> Material Properties: High Chemical Resistance

Strategic Partner of 

**CHEMTRIX**



EKasic<sup>®</sup> SiC plate: 2.5 yr, 180 °C with 50 % NaOH

✓ No sign of material corrosion



Together with exploiting the efficient thermal properties of EKasic<sup>®</sup> SiC, users also employ this material in harsh environments, for example;

- Nitrations
- Oxidations
- Chlorinations, Brominations & Fluorinations
- Wolff-Kishner reductions
- Alkylations



Suitable for control of exothermic processes

# EKasic<sup>®</sup> Material Properties: High Chemical Resistance

## Corrosion resistance of EKasic<sup>®</sup> Silicon Carbide



# Plantrix<sup>®</sup> Industrial Flow Reactor: Customer Application



## Reaction Challenges:

- Biphasic
  - Competing dinitration & decomposition products
  - Corrosive media
  - Challenging product isolation
- 
- Initially the reaction was investigated in a series of tube reactors (as illustrated)
    - A need for continuous mixing was identified Ekasic<sup>®</sup> SiC reactors selected




# Plantrix<sup>®</sup> Industrial Flow Reactor: Customer Application



Through the use of an intensified reactor, the reactant concentration & reaction temperature were increased – reducing the reaction time < 1 min

- Reduction in PMI (process mass index) across development & production



In addition to the reaction, controlled neutralisation & quenching is often key to maintaining reaction selectivity – generating waste – this can be reduced in flow

# Plantrix<sup>®</sup> Industrial Flow Reactor: Customer Application

**CHEMTRIX**



Strategic Partner of



DSM uses ESK Micro Reactors made of EKasic<sup>®</sup> (SiC) in a pharmaceutical production plant

# 'Coriac' Project

**Goal:** Benchmarking of multiple flow reactor technologies for a wide range of processes



Multipurpose pilot flow reactor skid for chemical conversion of suspensions

**Project Management:** TNO    **Timeframe:** End 2014

# 'Flow4API' Project

**Goal:** Explore 'Novel Process Windows' for API production in continuous flow

**Project Partners:**



**TNO**

Synthon

**CHEMTRIX**

*Cooperation between industrial partners and knowledge institutes*

**Focus:** Synthetic processes that pose a challenge for pharmaceutical, flavours & fragrances, veterinary and fine chemical industries. The project will focus on processes;


Involving High-value Molecules:

- Target yield improvements

Telescope Reactions:

- Cost savings for multi-step reactions

**Test Cases:** Provided by Synthon BV, explored using commercially available equipment, including the scalable flow platform from Chemtrix BV



**Project Management:** TNO (Gerhard Reeling Brouwer: [gerhard.reelingbrouwer@tno.nl](mailto:gerhard.reelingbrouwer@tno.nl))

**Timeframe:** End 2014

# Innovative Technology: Flow Reactor Benefits

**CHEMTRIX**

## 1. Safe Use of Extreme Reaction Conditions

- Efficient mixing
- Excellent thermal control
- Process intensification of hazardous reactions

## 2. Reduced Development Time

- Small hold-up volume
- Rapid reaction optimisation
- Minimal scale-up steps

## 3. Improved Process Control

- High level of reaction control
- Process reproducibility
- Quality by Design (QbD)

## 4. Reduced Production Costs

- Increased product quality
- Reduced safety investments
- Higher unit productivity

- ✓ Efficiency
- ✓ Quality
- ✓ Safety
- ✓ Sustainability

# Contact Details



Dr Charlotte Wiles (CEO)  
Chemtrix BV – Headquarters  
Chemelot Campus Gate 2  
Urmonderbaan 22  
6167 RD Geleen  
The Netherlands

e-mail: [c.wiles@chemtrix.com](mailto:c.wiles@chemtrix.com)

Tel: +31 4670 22600  
+44 1482 466459

Please find details of our publications, application notes and white papers  
on our website [www.chemtrix.com](http://www.chemtrix.com)