



micro
innova

efficient processing

Continuous Manufacturing:
Process Intensification Strategies
in Synthesis, Workup and Formulation
with a special focus on solids

Dr. Dirk Kirschneck, Microinnova Engineering GmbH



micro
innova

Content

- **Microinnova** Overview
- **Development** Strategy for **Flow** Processes
- Characteristics of **Flow Plant Design**
- Case Study: **Propoxylation**
- Case Study: **Crystallization**
- Case Study: **API-Plant**



efficient processing

www.microinnova.com



Microinnova Overview

7 Differentiators of Microinnova

1 | Chemist and Engineer Interaction

2 | Regulated Environment Experts for Continuous Manufacturing



3 | End-to-End Continuous Manufacturing Competence

- Synthesis
- Work-up
- Formulation



4 | Multi Process Intensification Technologies

4

6 | Tons per hour Competence



5 | 200+ Projects/WPs Experience



7 | Flow Processing of Solids



7



micro
innova

specialists in process intensification

**process
development**



**engineering
& plant**



**We focus on our customer's process
and provide high level solutions
independent from any
technology or supplier**


efficient processing

www.microinnova.com



micro
innova

MicroInnova's Approach to Success



The entire **process design** and **development work** until the final **turnkey plant** comprises of consecutive phases that enable a stepwise transfer of the existing process to a continuously realized process. A first plant sketch and cost estimation will be given after **Phase I**.

efficient processing

www.microinnova.com



micro
innova

Phase I – Process Design

A successful **business case study** for a specific product requires

- Theoretical data of the process to rate the overall process performance
- Estimations of investment costs

These informations will be available **after the Process Design Phase I.**

Theoretical Evaluation

Basic Lab Tests

Process Flow Diagram

Risk Assessment

Plant Cost Estimation

Phase II
Feasibility



efficient processing


www.microinnova.com

Focus on costs


Getting Costs down

Statement from BASF


1. Intensified Processes
2. Universal Engineering Designs



Better Control





New Strategies





Development Strategy for Flow Processes

How to design a flow process?





micro
innova

Drivers for Flow and/or Microreactors

mi




Yield

mi



Selectivity

mi



**Operation
Reduction**

mi



**Labor
Efficiency**

mi



Safety

mi



**Space-Time
Yield**

mi



Energy

mi



**Development
Speed**



efficient processing

www.microinnova.com



micro
innova

Application Profile




efficient processing

www.microinnova.com




micro
innova

MIC Fish Methodology



- k_1 ... desired standard reaction (solvent effect)
- k_2 ... follow up reaction (if late quench)
- k_3 ... side reaction (if $C_A \neq C_B$)
- k_4 ... **decomposition** (if high T)



- plug flow ?
- high T, p ?
- process liquid A ?
- process gaseous A ?




micro
innova

Technology Comparison




Engineer



efficient processing

www.microinnova.com

Process Intensification Toolbox



Micro
Mix



Micro Heat
Exchange




Unstable
Reagents



Plug
Flow



Cascade



Small
Volume



Mass
Transfer



Precise
Processing



Novel Process
Windows



Rotor
Stator




Quench



Microwave



Ultrasound



Membrane




Extrusion



micro
innova

Towards perfect processing




- narrow residence time distribution
- no backmixing
- no hot spots

- ideal stoichiometry
- no high concentration spots
- no dead zones

efficient processing

www.microinnova.com

QbD: Critical process parameter map



Out of the Box

Liquid-Liquid Processes

up to 400 °C

High Viscosity

Selfignition by O₂

High Corrosivity

Cryogenic Reactions

High exothermic reaction

Customer examples

5 out of Pharma TOP 10
 4 out of Generic TOP 10
 Chemical Global Players
 Polymer Companies
 Crop Protection

Gas-Catalytic Processes

up to 50 bar
 up to 300 °C

Liquid-Solid Processes

Precipitation
 Crystallization

Particle Surface
 Modification

Continuous
 Suspensions
 (cat/solid educt)

Melted
 Educt

Liquid-Gas Processes

Reactions
 with
 Cl₂/HCL

Reactions
 with O₂/H₂

Ozone
 (plant)

Reactions
 with NH₃ (liq)


**Very difficult
 Processes**



Process Design Strategies

Example: selectivity increase

Friedl-Crafts-Alkylation



driver



Selectivity



Yield

J. Yoshida, "Flash Chemsitry", Wiley, 2008

toolbox



Micro
Mix




Plug
Flow



Quench

Example: Unstable reagents



driver



Yield



A

k_1

B

k_2

C



toolbox



Precise Processing



Micro Mix



Plug Flow



Quench


Heejin Kim, Aiichiro Nagaki & Jun-ichi Yoshida




Nature Communications 2, Article number: 264




micro
innova

Value by optimizing mass & heat transfer



Operation	Batch	Flow	Factor
Reaction 	20 min	20 min	1
Mass Transfer 	50 min 1 m	0,5 min 1 cm	100
Heat Transfer 	25 min 0,5 m	0,5 min 1 cm	50



efficient processing



micro
innova

Technology Example: Extrusion

Continuous generation of one substream of a Healthcare formulation

Pharma Top10 Company



efficient processing

www.microinnova.com





micro
innova

Polymer Projects

- Free Radical Polymerisation
- Modification of Functional Groups
- Cross Linking
- Endcapping
- Oligomer Synthesis
- Formulation (e.g. Gel Structure)
- Encapsulation (Interfacial area Polymerization)



Engineering



- basic engineering
- detailed engineering
- 3D design
- automation solutions
- plant construction
- commissioning
- CE, ATEX, UL, UL-Ex, cGMP




Comparison batch versus conti / flow

advantage batch

- flexibility
- multipurpose


process is adjusted
to the plant



advantage conti/flow

- process performance
- safety
- easy automation

plant is adjusted
to the process



**concept necessary, which combines
batch flexibility with continuous
performance**




micro
innova

On-module flexibility

on-module adaption by exchanging specific parts

engineered spaces for adaption



efficient processing

www.microinnova.com

Modular plant design



utilities



automation + supply

- bus /plc
- electricity
- control air
- nitrogen


heating + chilling



www.huber-online.com




System Architecture Automation



Scale up of flow processes

Corning

AFR




Fluitec

XR



ESK

MR





IMM

Star-Lam



Flow Miniplant Level up to 10 l/hour




Manufacturing Level up to 10.000 l/hour



Case Study: **Propoxylation**


Case Study Propoxylation



Case study: Flow Miniplant

Flow Miniplant

Example of a 20 kg/h development or small scale production system





Case Study: **Crystallization**



micro
innova

Preferred particle size


Flame retardant

upper limit: mechanical properties of polymer
lower limit: retardant behaviour

Soluble salt

upper limit: speed of solubility
lower limit: dust formation

Cooling Crystallization



Case Study Flow Miniplant

Continuous Crystallization



Continuous Crystallization


Model Predictive Control



DN-15 reactor predictive controller


Multi-API Approach by Modules

Module 1 Module 2 Module 3 Module 4 Module 5




1. Starting Material

By changing starting material, different products with similar core structures can be obtained using the same modular sequence.




2. Reagent Choice

By changing the reagents within a module it is possible to access different classes of compounds.



3. Order of Modules

By changing the order of modules, different families of structures can be obtained





micro
innova

Continuous Cryogenic API Plant

Reaction: Continuous Solid Dosing into corrosive liquid under cryogenic conditions

- **Residence time improved:**
 - from 7 hours to 15 min
- **Process safety increased**
- **Lossless scale-up**



efficient processing

www.microinnova.com

Unit operation: Filter-dryer

Continuous Filtering & Drying







micro
innova

Services & Results

Service Map & Results



Results



Documents

- conceptual engineering
- basic engineering
- detailed engineering



-unit operation

Turnkey Plants




-modular/multipurpose



Reports

- consulting
- feasibility
- optimization



Lab experiments

- feasibility
- optimization




-redesigned


efficient processing

www.microinnova.com

Summary



**process
performance**



flexibility



safety



**fast track
realization**

**turning process
performance into money**

**efficient
processing**



mi

**micro
innova**

**process design
engineering
manufacturing plant**



Please contact me:

+43-3182-62626-110

dirk.kirschneck@microinnova.com