


*Enabling Catalysis in the Pharmaceutical Industry
with High Throughput Screening (HTS)*

Amazing where you can go

Asymmetric Hydrogenation: What It's All About



Metolachlor Josiphos Application (Syngenta's (S)-Dual)



Syngenta's Racemic Switch: Lifecycle Management of an Herbicide

- Originally marketed as a racemate (One enantiomer of the herbicide is active)
- 79% e.e. is OK for a herbicide
- 1,000,000:1 S/C ensures a cheap product with more potency than racemate
- 10,000 ton/year (10 ton batches, 16g Ir precursor & 32g ligand!)

Asymmetric Hydrogenation: Cheap, efficient, *surmountable* manufacturing technology


It's All About Options....



Choose the right lead for development!

- Productivity and activity (s/c and reaction time)
- Substrate quality tolerance
- Ligand/catalyst availability in bulk
- Price of catalyst (metal and ligand)
- Infrastructure requirements (alloy, pressure...)

Success Factors in Screening



HTS Equipment



Symyx core module in N₂ filled glovebox



96-well plate (0.4 mL, 100 bar H₂ or CO)

Highlights

- Inert handling
- Robotic dispensing
- SFC Analysis
- Symyx software


Benefits

- No O₂ (Piece of mind)
- Accuracy with robotics (Piece of mind)
- Fast time to hits
- Serendipity

Data Organization Made Easy: Asymmetric H₂ 96-Well Plate Analysis Example



LEA Software



The Staff



- All staff has an 'ausbildung' in automation
- All HTS - All day - Every day
- Consistency of operation (wheel not reinvented)
- No need to train every chemist on equipment
- Excellent suggestions to solve uncommon problems

The Ligand Library



- >600 Chiral Ligands
- (>100 CX Coupling ligands/catalysts)
- Storage (air sensitive/non-air sensitive)
- Updated on a regular basis

The Workflow

HTS plate design is done by a catalysis expert

→ Combination of know-how and (!) serendipity results in powerful shortcuts to solutions

Ligands, additives, metals: Test wide experimental space

→ Robot takes over dispensing and setting up of reactions
 < 10% of manpower and time required vs non-robotized setup

In situ Catalyst generation from ligand and precursor

→ Highest possible flexibility (metal, precursor type, ligand, counterion)

HTS Analytics

Need for speed (rt << 15 min / sample required)

→ SFC, GC, fast HPLC



The Know-How (Nuts and Bolts)




Applied research

- Comparison: hand exp's vs plate exp's.
 - How to stir (magnetic, shaking, stir shape)
 - Addition order
 - Residency times
- Cleaning robot needle after drawing a slurry.
- Dosing a solid.
- What's easy / what's hard?

Nuts and Bolts: The Solutions



Solids handling
(trap door unit)




Needle cleaning
(after drawing a slurry of hetero catalyst)



Filtration x 96

Case 1: Hydrogenation of Ketone



Targets:

- >90% ee
- s/c >1000


State of the Art:

- Only ca. 40% ee (PP)*(NN)*Ru Noyori-type catalysts
- Low e.e. likely due to 'catalyst poison' in R2 group

Design and Results of 1st HTS Plate (Rational and Serendipity Design)


- Precursors: $[\text{Rh}(\text{nbd})_2]\text{BF}_4$, $[(\text{Rh}(\text{nbd})\text{Cl})_2]$, $[\text{Ir}(\text{cod})\text{Cl}]_2$, $[\text{Ru}(\text{p-cymene})\text{I}_2]_2$
- Ligands: 16 ligands from 6 ligand classes
- 3 solvents and Et_3N additive

Best results



T001-1 (Taniaphos)

MeOH
 $[\text{RuI}_2(\text{p-cymene})]_2$
81% e.e.



A109-2 (MeO-Biphep)


MeOH/ Et_3N
 $[\text{Rh}(\text{nbd})\text{Cl}]_2$
87% e.e.

- A good start, but not quite good enough.
- 2nd plate will be designed based on results of the first emphasizing:
 - Taniaphos derivatives/Ru systems
 - MeO-Biphep derivatives/Rh systems


Design and Results of 2nd HTS Plate (Built Upon Results of the 1st HTS Plate)

- Precursors (all Ru and Rh): $[\text{Ru}(\text{p-cymene})\text{Cl}_2]_2$, $[\text{Ru}(\text{p-cymene})\text{I}_2]_2$, $[\text{Rh}(\text{nbd})\text{Cl}]_2$, $[\text{Rh}(\text{nbd})\text{TFA}_2]$, $[\text{Ru}(\text{cod})\text{TFA}_2]$
- Ligands: 11 ligands (incl. 8 from MeO-Biphep and Taniaphos classes)
- MeOH (single solvent) and Et_3N additive

Best results



T001-1 (Taniaphos)
MeOH/ Et_3N
 $[\text{Ru}]_2(\text{p-cymene})_2$
90% e.e.




T001-1 (Taniaphos)
MeOH/ Et_3N
 $[\text{Ru}(\text{cod})(\text{TFA})_2]$
85% e.e.

- Higher e.e.'s with the 'best' ligand from 1st plate (just different parameters!)
- Big advantage to screen ligands with multiple precursors of same metal (Ru)

Optimized to 3000:1 S/C with e.e.'s of 93%

Case 2: Hydrogenation of Acrylic Acid




Targets:


- >95% ee
- 3 kg product delivery ASAP

1st HTS plate:

- Highly diverse design of 32 ligands x 3 parameters (Rh, Ru precursors)
- Hits of 90% ee with Rh/P053-2 and Rh/F056-1:



kephos




fengphos

Design and Results of 2nd HTS Plate (Built Upon Results of the 1st HTS Plate)

- Precursors: $[\text{Rh}(\text{nbd})_2]\text{BF}_4$, $[(\text{Rh}(\text{cod})\text{Cl})_2]$, $[(\text{Ir}(\text{cod})\text{Cl})_2]$, $[\text{Ru}(\text{p-cymene})\text{I}_2]_2$
- Ligands: 27 ligands (incl. 16 from Kephos and Fengphos classes)
- 3 solvents: MeOH, DCE, Toluene

Best results




- Screening deeper into Fengphos and Kephos was unsuccessful
- BUT...Serendipity provided 2 hits (same ligand, different metals, different enantiomers)


Ru system optimized to 1000:1 S/C and scaled to 3 kg (5 wks!)

Tekturna: It all started with screening


Route design based on chiral synthons and known chemistry




'The CruX' is the development of an asymmetric hydrogenation



Solvias Catalysis Solutions: From Screening to Implementation




- Wide variety of commercial ligands screened
- Interesting push-pull ee trend shown in Walphos series w/ $[\text{Rh}(\text{NBD})_2]\text{BF}_4$ precursor:



- Screening + ability to systematically probe electronic/steric differences is powerful combination
- The advantage of a large, easily modifiable ligand library.

Industrial Chiral Ligands: Kg Scale Quantities for Your Process

Solvias Chiral Ligands




Josiphos


Walphos

Taniaphos
(Umicore)


Mandyphos
(Umicore)




MeO-Biphep



Butiphane




Naud-Cat.




Ubaphox


NEW from Evonik...




catASium D



catASium T



catASium MN



catASium M


Solvias has acquired an exclusive license from Evonik Degussa GmbH to manufacture, market, and use the **catASium** and **cataCXium** ligand product lines.

Synthon A key step for Tekturna process




Conditions:

- Walphos ligand
- $[\text{Rh}(\text{NBD})_2]\text{BF}_4$
- 100 bar H_2
- 80 °C
- 15 h
- **5000:1 S/C**
- **95% e.e.**



- Screening lead optimized for catalyst loading and manufacturing robustness
- Technical transfer of the best process to a 3rd party for manufacturing
- Process is performed at ton scale (with Solvias Walphos ligand)


CX Coupling: Widespread Industrial Prevalence



- Widespread industrial prevalence
- Tend to be developed with in-house expertise with standard reaction equipment
- Applying HTS screening alleviates 'lead bottleneck'


Focus on Industrial CN Coupling

Misc. Published Examples




24 Ligand x 4 Parameter Design (CN Coupling)


Diphosphines




Dppf
Generic




Dppp
Generic




Xantphos
Generic



Binap
Generic




Josiphos
Solvias




Josiphos
Solvias


Complexes



CC02
Solvias




CX31
Umicore




CX22
Umicore





Monodentate




cataCXium® A
Solvias




cataCXium® P series (Solvias)




PCy3
Generic




PPh3
Generic




P(o-anisyl)₃
Generic




P(MOD)₃
Generic




P(t-Bu)₃
Generic




Johnphos
Shasun




Sphos
Shasun



Xphos
Shasun




Ruphos
Shasun



c-bridp
Takasago

What is the Best Ligand System for CN Coupling?




- Buchwald system with X-Phos scaled to 100g at Saltigo (X-Phos, $\text{Pd}_2(\text{dba})_3$, NaOtBu, Toluene, MeOH, 78 °C).*

Can HTS provide alternative answers to this system?

*Presented by Dr. Buchwald at Dr. Pfaltz' 60th birthday symposium in 2008.


HTS Screening of Buchwald System




- Solvias performed standard 24 ligand x 4 design on transformation
- Judged conversion vs internal standard (because of sm volatility)
- Examined % conversion and by-product formation (incl. double addition)
- Parameters:
 - Pd(OAc)₂, Toluene, NaOtBu
 - Pd₂(dba)₃, DME, K₃PO₄
 - Pd(OAc)₂, Toluene, K₃PO₄
 - Pd₂(dba)₃, DME, NaOtBu
- Conditions: 1 eq. ArCl, 1.5 eq. Piperazine, 2 eq. Base, 0.25 eq. Trimethoxybenzene (internal standard), 80 C, 100:1 S/C, 16 hours

Screening results


Buchwald hits*



Pd₂(dba)₃, DME, K₃PO₄ (74%)
 Pd(OAc)₂, Tol, NaOtBu (66%)
 Pd₂(dba)₃, DME, NaOtBu (57%)



Pd₂(dba)₃, DME, K₃PO₄ (67%)
 Pd(OAc)₂, Tol, K₃PO₄ (56%)




Pd₂(dba)₃, DME, K₃PO₄ (69%)


•<0.5% double addition observed for these hits!

•Highest conversion for generic ligand was 47% (system: t(Bu)₃P)


Other hits*




DME, K₃PO₄ (63%)



Pd(OAc)₂, Tol, NaOtBu (72%)
 Pd₂(dba)₃, DME, K₃PO₄ (61%)




Pd₂(dba)₃, DME, K₃PO₄ (70%)



Pd(OAc)₂, Tol, NaOtBu (57%)

HTS Buchwald substrate conclusion



Four alternative viable ligand systems were found!

-Allows for the choice of the *best* system for development (chemistry and cost)

Acknowledgements



- Heidi Landert
- Hans Meier
- Dr. Xiangdong Feng
- Dr. Felix Spindler
- Dr. Ulrike Nettekoven
- Dr. Martin Kesselgruber
- Dr. Richard Sachleben

- Dr. Fred Naud
- Dr. Matthias Lotz
- Dr. Pierre Martin
- Dr. Marc Thommen
- Dr. Anita Schnyder
- Dr. Willy Lottenbach
- Jeannine Müller
- Markus Müller
- Heinz Steiner
- Urs Siegrist
- Reto Kellerhals

- Dr. Hans-Jochen Dannappel



Prof. Walter Weissensteiner, University of Vienna
Prof. Armin Börner, Ifok Rostock
Prof. Albrecht Salzer, RWTH Aachen
Prof. Olivier Riant, Université Louvain La Neuve

Prof. Andreas Pfaltz, University Basel
Prof. Antonio Togni, ETH Zürich
Prof. Pher Andersson, University of Uppsala
Federal Commission for Technology & Innovation